



Demonstration of blood A and Rh blood group of a blood sampling held on September 1

Aim:

To determine the ABO and Rh blood group of a blood sample using the slide/tile method.

Principle

Blood grouping is based on **agglutination (clumping)** of red blood cells when they react with specific antisera

Anti-A → detects A antigen

Anti-B → detects B antigen

Anti-D → detects Rh (D) antigen

Note: This is a laboratory demonstration. For clinical transfusion or medical purposes, blood grouping should be performed using validated laboratory methods and appropriate controls.

